

porodo



(Porodo Smart Wearable Ring Size 8 &10
SKU: PD-SRG10 **SKU: PD-SRG8**)

Table of Content

App User Manual_____	2
Installing the ChipletRing App_____	2
1. Android App_____	3
1.1 Personal Information_____	3
1.2 Device Connection_____	4
1.3 Measuring Heart Rate and Blood Oxygen__	5
1.4 Activity Monitoring_____	6
1.5 Sleep Monitoring_____	7
1.6 Additional Features_____	8
2. IOS App_____	10
2.1 Personal Information_____	10
2.2 Device Connection_____	11
2.3 Measuring Heart Rate and Blood Oxygen__	12
2.4 Activity Monitoring_____	13
2.5 Sleep Monitoring_____	14
2.6 Additional Features_____	15
Warranty_____	17
Contact Us _____	17

App User Manual

Installing the ChipletRing App

This document provides information about the Porodo Smart Wearable Ring. Before starting to use the product, please read this manual carefully. Scan the following QR codes to download and install the "ChipletRing". Upon downloading and installing the app, enable Bluetooth connection on the smartphone. Search for the device name "Por-Ring". Connect to the device.



iOS




Android

1. Android App

1.1 Personal Information

Start by downloading and installation of the ChipletRing App. Next, sign up by filling out the required personal information.





The screenshot shows the 'Profile' screen of the ChipletRing app. The status bar at the top displays the time 17:15, signal strength, Wi-Fi, and battery icons. The app bar has a back arrow and the title 'Profile'. Below the app bar, there are three input fields for personal information: Height (175 cm), Weight (65 kg), and Step (80 cm). Each field consists of a label, a text input box, and a unit. The bottom of the screen shows a white home indicator bar.

Field	Value	Unit
Height	175	cm
Weight	65	kg
Step	80	cm

1.2 Device Connection

To add the device, enter the App's home page and choose "Add Device."



If several devices are recognized, check the MAC ID of the product and choose from the list accordingly.



1.3 Measuring Heart Rate and Blood Oxygen

Select the marked items as shown in following figures to measure

1. Heart rate
2. Blood Oxygen saturation
3. Heart rate fluctuations
4. Body temperature





1.4 Activity Monitoring

Every time the user uses the device, their activity is logged, and the data may be viewed for the specified day, week, or month.




1.5 Sleep Monitoring

To view the data recorded during sleep, select the marked item. Note that the data stored in App is updated each time the Smart Wearable Ring is connected to ChipletRing App.



1.6 Additional Features

To set specific activity goals and change the basic information, select the small icon indicated in the figure.



To set activity goals change the settings in the App.

Note:

The user can change and update the basic physical information, including weight and height.


The screenshot displays the 'Goal' settings interface of a fitness application. At the top, the status bar shows the time as 17:18 and icons for cellular signal, Wi-Fi, and battery. Below the status bar, a dark blue header contains a back arrow on the left and the title 'Goal' in the center. The main content area is divided into three sections, each with a label on the left and a value input field on the right. The first section is labeled 'Steps' and has a value of '3500' with the unit 'Step' to its right. The second section is labeled 'Calories' and has a value of '3000' with the unit 'kCal' to its right. The third section is labeled 'Workout Time' and has two input fields: '0' for 'Hour' and '30' for 'min'. The bottom of the screen shows a dark blue home indicator bar.

Category	Value	Unit
Steps	3500	Step
Calories	3000	kCal
Workout Time	0	Hour
Workout Time	30	min

2. IOS App

2.1 Personal Information

Download and install the ChipletRing App. Open the App and fill out the required information.





The screenshot shows the 'Profile' screen of the ChipletRing App. The status bar at the top displays the time 17:15, signal strength, Wi-Fi, and battery icons. The screen has a dark blue background. A back arrow is in the top left, and the title 'Profile' is in the top center. Below the title, there are three input fields for personal information, each with a label on the left and a unit on the right. The first field is for 'Height' with a value of '175' and unit 'cm'. The second field is for 'Weight' with a value of '65' and unit 'kg'. The third field is for 'Step' with a value of '80' and unit 'cm'. The bottom of the screen shows a white home indicator bar.

Field	Value	Unit
Height	175	cm
Weight	65	kg
Step	80	cm

2.2 Device Connection



To add the device, enter the App's home page and choose "Add Device."

If several devices are recognized, check the MAC ID of the product and choose from the list accordingly.



2.3 Measuring Heart Rate and Blood Oxygen

Select the marked icons to measure the heart rate, blood Oxygen saturation, and changes in the heart rate.





2.4 Activity Monitoring

To check the activity log, select the "Activity" item. The user can observe the activity data for a specific day or desired weeks and months.




2.5 Sleep Monitoring

This device is capable of automatically recording data during the time the user is sleeping. Select the "Sleep" icon to preview the sleep data.



2.6 Additional Features

Select the "Me" item to set activity goals or to change personal information.



To set activity goals change the settings in the App.

Note:

The user can change and update the basic physical information, including weight and height.

The screenshot displays the 'Goal' settings interface of a fitness application. At the top, the status bar shows the time as 17:18 and various system icons. Below the status bar, a dark blue header contains a back arrow and the title 'Goal'. The main content area features three distinct settings sections, each with a label on the left and a white input field on the right. The first section is for 'Steps', with the value '3500' and the unit 'Step'. The second section is for 'Calories', with the value '3000' and the unit 'kCal'. The third section is for 'Workout Time', which includes two input fields: '0' for 'Hour' and '30' for 'min'. The bottom portion of the screen is a solid dark blue area, and a white home indicator bar is visible at the very bottom.

Category	Value	Unit
Steps	3500	Step
Calories	3000	kCal
Workout Time	0	Hour
Workout Time	30	min

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **<https://porodo.net/warranty>** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:

<https://www.porodo.net/warranty>

Contact Us

If you have any questions about this Privacy Policy, please contact us at: **info@porodo.net**

Website: **porodo.net**

Service Support: **support@porodo.net**

Instagram: **[porodo](https://www.instagram.com/porodo)**