

PORODO
LIFESTYLE



Porodo Lifestyle
SMART SONIC
ELECTRIC TOOTHBRUSH

SKU: PDLFSTXCA8



Table of Contents

Product Overview	2
Specifications	2
Setup Guide	3
Disposal	5
Warranty	5
Contact Us	5

Product Overview

The Porodo Smart Sonic Electric Toothbrush offers a powerful and customizable cleaning experience with six cleaning modes to suit different needs, including Clean, Mild, Whitening, Polish, GumCare, and Sensitive. With IPX7 water resistance, it is safe to use in the shower, and the toothbrush provides up to 30 days of use on a full charge. The built-in smart timer and reminder ensure an efficient brushing routine, while the Dupont bristles guarantee a high-quality cleaning. The toothbrush also features a memory function for convenience.

Specifications

Input	5V/1A
Battery Capacity	500mAh
Charging Time	4h
Working Time	Up to 30 Days
Water Resistance Level	IPX7
Charging Interface	USB-C
Product Weight	64.5g
Product Size	30 × 30 × 233mm
Model Number	PDFLSTXCA8

Setup Guide

A. Instructions:

- 1.** Squeeze an appropriate amount of toothpaste onto the brush head and wet it (ensure the toothbrush is not turned on).
- 2.** Gently press the toothbrush head against your teeth.
- 3.** Open the toothbrush cover (Note: open it after placing it in your mouth).
- 4.** Move the toothbrush slowly, cleaning all areas of the teeth as per the tips.
- 5.** After the toothbrush stops, remove it and rinse both the toothbrush body and brush head with water.

B. Brush Head Description:

- 1.** Installation: Align the bristles with the button, insert them directly with slight force, and avoid twisting the brush head to prevent the main body from breaking.
- 2.** Disassembly: If the bristles are aligned with the button, pull them out directly with a little force, ensuring you do not twist the brush head, which may cause the main body to break.
- 3.** If there are gaps between the brush head and the main body of the toothbrush, this is normal and does not affect operation.
- 4.** To achieve better whitening results, it is recommended to replace the brush head once a month.
- 5.** If the bristles begin to bend outward, it is time to replace the brush head.
- 6.** Always use the recommended brush head (please purchase from the store where you bought the toothbrush).
- 7.** The bristles may deform during production or transportation. This is normal and can be resolved by soaking them in warm water for 10 minutes.

C. Toothbrush Charging Instructions:

- 1.** Remove the waterproof sleeve and insert the DC round interface charging cable.
- 2.** Plug the charging cable into a suitable power supply (e.g., computer, power bank, USB charging plug).

3. The toothbrush uses a low current design, which can be continuously charged for more than six hours.
 4. Pull out the charger after full charge.
 5. Avoid charging in high-temperature environments, as it may cause overheating or component damage.
- D. Matters Needing Attention:**
1. After using the electric toothbrush, you may find slight bleeding in your gums. This phenomenon is normal and usually disappears in a few days. If your gums continue to bleed, consult your dentist.
 2. If you have had surgery on your mouth or gums in the past two months, please consult your dentist before using this toothbrush.
 3. This electric toothbrush is not suitable for the following groups (including children under 6 years of age): persons with physical function or intellectual degradation, as well as those who lack experience and use knowledge. These people can only use the electric toothbrush after having a person to ensure the personal safety of others nearby to supervise or guide the use of the electric toothbrush.
 4. Do not use hair care products (such as shampoo, shower gel, hand washing liquid, gel, oil) to contact the product so as not to damage the surface of its product.
 5. Do not use a toothbrush head other than that recommended by the manufacturer (please purchase the toothbrush head in the store where the product has been purchased).
 6. Applicable temperature of toothbrush host ($0^{\circ}\text{C} < 40^{\circ}\text{C}$); Applicable temperature of brush head ($0^{\circ}\text{C} < 80^{\circ}\text{C}$)
 7. This product is for private use only. Dental clinics or multiple people cannot be used together.
 8. This product is only designed to clean teeth and gums. Do not use for any other purpose.
 9. The toothbrush shall not be exposed to the direct sun, but shall be placed in a dry and ventilated place at $0^{\circ}\text{C} < 40^{\circ}\text{C}$.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **porodo.net/warranty** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
porodo.net/warranty

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: **porodo.net**

Service Support: **support@porodo.net**

Instagram: **[porodo](https://www.instagram.com/porodo)**