

PORODO
LIFESTYLE



Porodo Lifestyle

**MINIMEASURE
ELECTRIC DIGITAL BABY SCALE**

SKU: PDLFSTDC5HWH

Table of Contents

Product Overview	2
Specifications	2
Package Contents	3
Buttons and Functions	3
Installing the Batteries	5
Usage Instructions	5
Safety Precautions	6
FCC Statement	7
Disposal	8
Warranty	8
Contact Us	8

Product Overview

The Porodo Lifestyle Minimeasure Electric Digital Baby Scale is a precise and reliable weighing scale designed specifically for infants. It features a high-precision sensor that ensures accurate measurements up to 40kg, making it suitable for babies and young children. The scale includes a 150cm measuring tape for easy height measurement and offers unit conversion capabilities for added convenience. The tare function allows for zeroing out the weight of blankets or other materials, while the hold function ensures the measurement remains displayed for easy reading. This, user-friendly scale is ideal for parents who look for a trusted way to monitor their child's growth.

Specifications

Material	ABS
Units	g, kg, oz, lb
Power Supply	3xAAA Batteries
Weighing Range	30g-40kg
Division	5g
Length Measuring	0cm-150cm
Product Weight	1.05kg



Product Size	47 x 320 x 550mm
Model Number	PDLFSTDC5H

Package Contents

Scale x1

AAA Batteries x3

Buttons and Functions



1. ON/OFF Button:

Press this button to turn the scale on, and press again to turn it off.

2. "Tare" Button:




- a. Press the "TARE" button to clear the scale's own weight before weighing.
- b. Press "TARE" before placing clothes and blankets on the tray to remove the weight of the clothes.
- c. Once the display shows "0" (as shown below), place the baby on the scale to ensure an accurate measurement.

3. "Hold" Button: Hold Weight

- a. Once the weight value has become stable, press the "HOLD" button to capture the weight on the screen.
- b. To cancel the captured data, press the button again.
- c. The weight will automatically lock. Once fixed, the weight will flash twice and then lock.
- d. To hold the weight for an extended period, press the "HOLD" button again.


4. "UNIT" Button: Select Unit

Use the "UNIT" button to select the unit of weight displayed on the screen. You can measure in different units like kilograms (kg) or ounces (oz).



5. Automatic Power Off Function (Energy Saving Mode):

The scale automatically switches off after approximately 2 minutes of inactivity to save battery.



6. Battery Consumption:

This product is energy-efficient and environmentally friendly. It operates with 3 AAA batteries.

Ensure the insulation film is removed from the batteries before use, and install them with the correct "+" and "-" orientation.

Installing the Batteries

1. Open the battery cover located at the bottom of the scale.
2. When installing the batteries, place the "-" side of the batteries into the battery compartment first, then press the "+" side into place to complete the installation.

Note: Ensure the correct orientation of the "+" and "-" sides of the batteries before installation.



Usage Instructions

1. Place the scale on a table, desk, or any other flat surface. It should not be placed on a carpet.



2. Press the "ON/OFF" button to turn the scale on, then press the "TARE" button to ensure the display shows "0".
3. Once the display shows "0", you can start using the scale.



4. After the display shows "0" (as shown in the picture), place the baby on the scale. Press the "UNIT" button to select your preferred unit (g, kg, oz, lb: oz).
5. Place the baby on the weighing plate with its full weight. The digital display will show the weight in the unit you have selected.
6. Press the "HOLD" button to freeze the weight display. Press the button again to delete the retained data.

NOTE: If you are using the scale without the weighing plate, keep the baby upright on the scale and ensure the baby's feet are evenly placed so the weight is distributed properly.

Safety Precautions

1. Place the scale on a flat surface before use.
2. Do not place or use the product in a damp, very hot, or very cold environment.
3. Protect the product from being damaged by falls or strong shocks, as it is a high-precision electronic product.
4. Only use a soft cleaning cloth; do not use chemical solutions, water, or other liquids to clean the product.
5. This product is intended for measuring family data and should not be used for professional measurements.
6. Remove the batteries if you are not going to use the product for a long period.
7. Replace the batteries when the LED display shows "L0".

8. If the LED display shows "ERR" (indicating overload, with a maximum weight capacity of 40kg), please remove any objects to avoid damage.

FCC Statement

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference; and

(2) This device must accept any interference received, including interference that may cause undesired operation.

2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

Reorient or relocate the receiving antenna.

Increase the separation between the equipment and receiver.

Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

Consult the dealer or an experienced radio/TV technician for assistance.

RF warning statement:

The device has been evaluated to meet general RF exposure requirements.

The device can be used in portable exposure conditions without restriction.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at <https://www.porodo.net/warranty> and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
<https://www.porodo.net/warranty>

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: <https://www.porodo.net/>

Service Support: **support@porodo.net**

Instagram: **porodo**