

PORODO
LIFESTYLE



Porodo Lifestyle

**Relieva
Roller Massager**

SKU: PDLFSTFH24BK

Table of Contents

Product Overview	2
Specifications	2
Product Components	3
Instructions for Use	3
Recommended Usage	4
Charging Guidelines	4
Safety Instructions	5
Disposal	5
Warranty	6
Contact Us	6

Product Overview

The RELIEVA Roller Massager is designed for a comfortable and relaxing massage experience. It features an ergonomic handle for easy use and six free-spinning rollers that provide a soothing massage. With three adjustable speed levels, you can customize your massage to your liking. The device also has an auto-off function for safety, turning off after a set time. It offers up to 60 minutes of working time, making it perfect for extended use. This massager is ideal for those looking to relieve tension and enjoy a relaxing massage at home.

Specifications

Material	ABS
Input Voltage	7.4V
Power	35W
Rated power	10W
Battery Capacity	1800mAh
Speed Levels	3
Working Time	60 Minutes
Charging Time	3h
Charging Port	USB-C
Product Weight	815g
Product Size	75 × 168 × 200mm
Model Number	PDLFSTFH24

Product Components



A. Power/Speed Adjust Button

B. Power/Speed Indicator Lights

C. USB-C Charging Port

D. Multi-grip Handle

E. Roll Cage

Instructions for Use



1. Power On: Press and hold the power button on the top of the handle for about 1 second. The roll cage will start spinning automatically.

2. Set Speed: Press the power button briefly to set the speed level (1-3), as indicated by the lights.

3. Power Off: Press and hold the power button on the top of the handle until the device turns off.

Recommended Usage



Shoulders

Neck

Arms

Legs



Chest

Charging Guidelines

1. To charge your RPM Mini, plug the included USB-C charging cable into the charging port on the outside of the handle.
2. The indicator lights will show the level of charge.

Note: Wall adapter for the USB-C cable is not included.



USB-C Charging Port

Safety Instructions

1. The RPM Mini is for home use only and not intended for children.
2. Do not use the RPM Mini on your head, face, or near your genitals.
3. Consult your doctor before using the RPM Mini if you are pregnant.
4. Keep loose hair or clothing away from any moving parts of the RPM Mini.
5. Do not use the RPM Mini above your throat or near your C4 vertebrae.
6. Do not get the RPM Mini wet. Clean the device with a lightly damp cloth only.
7. Do not cover the device or block the motor vents.
8. Never attempt to disassemble the RPM Mini.
9. Never store the RPM Mini at temperatures above 40°C/104°F. Only charge the RPM Mini when the ambient temperature is between 0°C/32°F and 40°C/104°F. If the device feels hot during use, stop immediately.
10. Never use or charge the RPM Mini near water or in areas with moisture.
11. For long-term storage, ensure the RPM Mini is fully charged but not plugged into the charger.
12. This product is not intended to diagnose, treat, or prevent any diseases. If you experience discomfort or pain during use, immediately stop and consult a doctor.
13. While the RPM Mini is unlikely to cause injury, there is a small risk that it could worsen existing conditions or injuries. You assume all risks associated with the use of the RPM Mini.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.

Warranty

Products that you buy directly from our **Porodo** website or shop come with a 24-month warranty.

When you buy **Porodo** products from any of our approved sellers, you only get a 12-month warranty. If you want to extend this warranty, go to our website at **<https://www.porodo.net/warranty>** and fill out the form with your information. Don't forget to upload a picture of the product too. After we've checked and accepted your request, we'll send you an email to confirm that your product's warranty has been extended.

For more info, please check:
<https://www.porodo.net/warranty>

Contact Us

If you have any questions about this Privacy Policy, please contact us at:

info@porodo.net

Website: **<https://www.porodo.net/>**

Service Support: **support@porodo.net**

Instagram: **porodo**